

Proudly Canadian!

FEATURES

- Tailored docking station with centralized charging, available in 2, 4 and 6 modules.
- Recognize 5 ml, 10 ml, 20/30 ml, and 50/60 ml syringes from major brands automatically (open system).
- Wide range of flow rate from 0.1 to 1000 ml/h.
- Flexible infusion modes: ml/h, time/volume and weigh mode.
- Enhance safety by using drug library, reduce the risk of over dose, memory can store up to 120 drug dosing protocols.
- Full alarm protection.
- Lithium battery powers the pump up to 24 hours.



Rotate 90 degrees for different drug treatments





PRODUCT SPECIFICATIONS

Infusion Rate Range	5 ml syringe
Flow Rate Increments	0.1-100 ml/h at 0.1 ml/h increments 101-1000 ml/h at 1ml/h increments
Infusion Modes	ml/h volume/time drug library patient weight
Infusion Time	1 min-99 hours
Bolus Volume	Programmable from 1 to 50 ml.
Bolus Rate	Programmable from 0.1 to 1000 ml/h based on syringe size
Accuracy	± 5% (± 2% after syringe calibration)
KVO (keep vein open)	Automatically switches to KVO rate (from 0.5 to 1 ml/h) at the end of infusion process
Occlusion Pressure	5.8 to 14.5 psi (0.04 to 0.10 Mpa), adjustable from High, Medium, Low
Alarm Volume	Selectable from High, Medium & Low
Alarms	Syringe not Loaded, Start Reminder, Near Empty, Complete, Occlusion, Low Battery, Syringe Disengaged, Mechanical Error, Infusion Error, Battery Depleted
Event Log Capacity	Maximum 800 records of operational activities or alarms with date and time of the event
Power Supply	AC 100-240V 50/60 Hz
Battery	Built-in rechargeable Lithium battery Pump will operate up to 24 hours at 5 ml/h with fully charged battery
Classification	Class IIb, Type CF, IPX2
Dimensions	21O(L) x 10O(W) x 45(H) mm
Weight	500g (battery included)



Canafusion Technologies Inc.

302-2999 Underhill Ave. Burnaby, BC Canada. V5A 3C2 Tel: +1(604) 568.0936 / +1(604) 568.6977 Fax: +1(604) 677.5295

canafusion.com • info@canafusion.com



