

### Delivering Nasal High Flow at home and in long-term care












The myAIRVO 2 provides a compact and effective solution for delivering the benefits of Nasal High Flow.





# Introducing the F&P myAIRVO 2, delivering respiratory support and airway hydration with a focus on comfort and ease of use.

Fisher & Paykel Healthcare's advanced humidification technology allows the myAIRVO 2 to comfortably provide high flows of air/oxygen mixtures to spontaneously breathing patients, through the unique Optiflow™ nasal cannula interface.





#### 1. Optiflow™ Nasal Interface

- · Soft, flexible prongs
- Wide-bore design reduces gas jetting
- Prongs contoured to the patient's nose

## 2. Supplementary oxygen (when required)

- Oxygen can be added from a concentrator or cylinder
- Inbuilt ultrasonic oxygen analyzer requires no calibration, service or replacement

## 3. Designed for simple setup, use and cleaning

- Helpful onscreen animations assist with setup and troubleshooting
- Quiet operation day and night

### 4. Adjustable temperature and flow settings

- Three temperature settings (37, 34, 31 °C) help achieve comfort and compliance
- Integrated flow generator delivers a wide flow range (up to 60 L/min) - no compressor or wall air supply required

#### 5. Heated breathing tube

- Dual spiral heater wires and unique integrated temperature sensor
- No separate temperature probes or heater-wire adapters required

### 6. An integrated system

 The F&P myAIRVO 2, with its wide range of accessories, provides versatility, mobility and convenience

